

# ProboStat™ ACIS

## Sample Holder for AC Conductivity Measurements and Impedance Spectroscopy

**ProboStat™ ACIS is a sample holder dedicated for AC conductivity measurements and 2-electrode impedance spectroscopy on small disk samples at high temperatures and under controlled atmospheres.**

Special design of the hot-zone parts facilitates mounting and replacement of test samples compared with the standard ProboStat™ furnishment.

The test sample rests in a support tube, is contacted by electrode contacts, and held stably in place with a spring load assembly. The setup is assembled in a single chamber mode and can be fed with virtually any gas. The sample holder is designed to operate at near-atmospheric pressure, but can be pumped to low vacuum, or hold up to 15 bar. Electrical connections are made via standard coax cables.

### ProboStat™ ACIS Compatibility

The system is supplied with a standard ProboStat™ base unit and can be expanded for any sample geometries and measurements that ProboStat™ can support.

The ProboStat™ ACIS works well with all instrumentation such as impedance spectrometers, LCR meters and potentiostats from Solartron, PAR, Novocontrol, Autolab, Gamry, HP/Agilent, etc.

### ProboStat™ ACIS Specifications

#### Temperature range:

- RT – 1600°C

#### Atmospheres:

- Oxidizing, inert, reducing
- Wet, dry
- Corrosive optional
- 0.01 mbar – 1 bar (15 bar optional)

#### Base unit:

- Ni-plated brass (SS316 optional)
- 6 BNC electrode connections
- 3 thermocouple connections
- 4 Swagelok quick-connects - two gas in and two gas out
- 16 mini-contacts for detachable electrodes and thermocouples
- 3 switches for grounding and shielding options
- Cooling/heating water hose fittings

#### Hot zone parts:

- Specially designed sample support tube of alumina
- Spring-load assembly of alumina
- 2 electrode contacts: Pt+alumina
- Thermocouple: type S or K

#### Test sample size:


- Disk diameter: max 15 mm
- Height: max 10 mm

#### Outer system dimensions:

- Outer tube diameter: 40 mm
- Overall length: 75 cm

### Materials properties measured and applicable methods

- AC conductivity vs T, pO<sub>2</sub>, pH<sub>2</sub>O, etc.
- Impedance spectroscopy
- Bulk, grain boundary and electrode impedance
- H/D isotope effect
- Fuel cell components test
- Dielectric properties



**Fig.1 ProboStat™ ACIS:**  
(a) hot-temperature zone  
(b) fully assembled

#### Hot-temperature zone:

1. Test sample
2. Spring-loaded system
3. Electrode connects
4. Sample support tube. The standard thermocouple is inside of the sample support tube

# NORECS products

## Measurement Systems

<b>ProboStat™</b>	Versatile sample holder system for measurements of electrical properties, transport parameters, and kinetic of materials, solid/gas interfaces and electrodes at high temperatures under controlled atmospheres.
<b>SeebSys</b>	System for automated measurements of Seebeck coefficient and electrical resistance at high temperatures and under controlled atmospheres. Measurement software included.
<b>CoultCell</b>	Coulometric titration system for measurements of materials oxygen nonstoichiometry at high temperatures under different oxygen partial pressure. Measurement software included.
<b>Tubular membranes and electrolytes for ProboStat™</b>	Use standard ProboStat™ for tubular gas separation membranes or electrolyte cells. Catalytic Membrane Reactor (CMR) for tubular ceramic membranes with CMR-modified ProboStat™.

## Software

<b>Omega</b>	Multipurpose measurement and data analysis software for electrochemical measurements. Omega is designed for maximum flexibility: it accepts data from many different sources such as instruments, furnaces and mass flow controllers. It can plot any measured or calculated data. With no time restrictions nor limit for simultaneous measurements. It was originally developed for ProboStat™ users.
<b>Omega Temp</b>	Software to control and program furnace temperatures and profiles, with unlimited number of segments.

## Atmosphere control

<b>ProGasMix</b>	Versatile rotameter-based manual gas mixer especially developed for laboratory investigations of the properties of oxides at high temperatures vs T, pO <sub>2</sub> , pH <sub>2</sub> O etc., or fuel tests with humidification of fuel and/or oxidant. Partial pressures are calculated via accompany software. Allows measurements of proton and oxide ion transport numbers separately.
<b>FCMix</b>	Simple and low-cost gas mixer that controls the flows of fuel, oxidant and one inert gas that can be routed to the fuel and oxidant for flushing, soft start, slow SOFC anode reduction, and safe operation, as well as tests of gas diffusion limitations.
<b>Probble</b>	A combined humidification unit and over-pressure controller suitable for small gas flows to SOFC button cell testing, annealing furnaces, etc. Two units used together with FCMix give combined humidification and overpressure control of each ProboStat™ cell chamber as well as fine pressure control in and between two cell chambers.

## We supply

<b>Furnaces</b>	Furnaces customized for ProboStat™
<b>TempBox</b>	Simple thermocouple reader with USB port
<b>PLD targets</b>	8YSZ, BZY-Ni, BZY, etc.
<b>SOFC button cells</b>	Anode supported cells with 8YSZ electrolyte, porous nickel cermet anode and lanthanum cobaltite cathode